	biokit HSV-2 Rapid	HerpeSelect®		CAPTIA™ HSV IgG Type
	Test	ELISA	Immunoblot	Specific ELISA
Supplier	biokit USA (formerly Diagnology with POCkit® test)	Focus Diagnostics, Inc. (formerly Focus Technologies)		Trinity Biotech USA
FDA Approved	Yes (August 1999)	Yes (February 2000)	Yes (April 2000)	Yes (July 2004)
Antibodies Detected	HSV-2 only	HSV-1 <u>or</u> HSV-2	HSV-1 and HSV-2	HSV-1 <u>or</u> HSV-2
Sensitivity (no false negatives)	93-96% accurate*	96-100% accurate*	97-100% accurate*	90-100% accurate**
Specificity (no false positives)	95-98% accurate*	97-100% accurate*	98% accurate*	90-99% accurate**
HRC Recommended Window Period	12 to 16 weeks after exposure	12 to 16 weeks after exposure	12 to 16 weeks after exposure	12 to 16 weeks after exposure
Collection Method	Finger prick/blood	Blood draw (sent to local lab)	Blood draw (sent to local lab)	Blood draw (sent to local lab)
Test Location	In provider office	Various labs	Various labs	Various labs
Result Time	Less than 10 minutes	~ 1 to 2 weeks	~ 1 to 2 weeks	~ 1 to 2 weeks
Can Be Used During Pregnancy	Not FDA approved	Yes	Yes	Yes
Manufacturer's List Price Per Test (does not include provider or lab fees)	~ \$20	~ \$4	~ \$19	~ \$4
Phone Number	Toll-free (800) 926-3353	Toll-free (800) 505-0536 (to find a specific lab)		Toll-free (800) 325-3424 (to find a specific lab)
Web site	http://www.biokit.com/	http://www.herpeselect.com/		http://www.trinitybiotech.com

Table 5. Herpes Blood Tests Quick Reference Guide (from http://www.ashastd.org/pdfs/blood_test.pdf, accessed February 3, 2005)

* Adapted from Reference (8).

** Adapted from Reference (8), and from FDA Summary of Safety and Effectiveness Information CAPTIA[™] HSV-2 IgG Type Specific Test Kit, July 2004.

Also available from Fisher HealthCare as Sure-Vue[™] HSV-2 Rapid Test. For more information, call toll-free (800) 766-7000 or log on to www.fisherhealthcare.com.